

A+ Aloe™ Nitrile Gloves

Formulated with skin-protecting aloe

Powder-free gloves feature a unique aloe coating on the interior glove surface for a comfortable, soothing feel against bare hands. Aloe prevents the drying and chapping that plagues frequent glove wearers.

Gloves are latex-free and contain no natural rubber proteins. Synthetic formulation offers good chemical and puncture resistance. Beaded cuff ensures easy donning. Ambidextrous.



86480-10

Specifications & Ordering Information

Catalog number	Size	Length	Thickness	Qty/bx	Price per box	
					1 bx	10 bx
K-86480-10	Small	9½"	5 mil	100		
K-86480-12	Medium					
K-86480-14	Large					
K-86480-16	X-large					

[K-86402-00](#) Wire glove dispenser, single

GO to page(s) 1718 to 1727

We offer a complete supply of lab coats, coveralls, aprons, sleeves, and boot covers. See our "Safety-Apparel" section.

MICR OFLEX® FreeForm™ SE Nitrile Gloves

Patented "super soft" formulation feels like latex

Powder-free, 100% nitrile gloves offer superior comfort and elasticity in a synthetic glove. Unique formulation conforms to the contours of your hands as you move—minimizing unnecessary hand stress and enabling a full range-of-motion.

Lightly textured fingers ensure a secure grip in wet and dry conditions. Interior polymer coating assists in quick and easy donning.

Compliance: ASTM D 6319.



86445-01

Specifications & Ordering Information

Catalog number	Size	Length	Thickness	Qty/bx	Price per box	
					1 bx	10 bx
K-86445-00	Small	9½"	6 mil	100		
K-86445-01	Medium					
K-86445-02	Large					
K-86445-03	X-large					

Ansell Protective Products Touch N Tuff™ Nitrile Gloves

The standard in comfort, durability, and cleanliness

Premium quality gloves are thin enough to allow excellent dexterity, yet tough enough to resist punctures and tears. Gloves are 100% nitrile and contain no fillers, silicones, or plasticizers. Ambidextrous.

Approval: CE

Compliance: FDA regulations for food contact



81602-55

Specifications & Ordering Information


Catalog number	Size	Length	Thickness	Qty/bx	Price per box	
					1 bx	10 bx
Powdered gloves						
K-81602-40	Small	9½"	4 mil	100		
K-81602-50	Medium					
K-81602-60	Large					
K-81602-70	X-large					
Powder-free gloves						
K-81602-45	Small	9½"	4 mil	100		
K-81602-55	Medium					
K-81602-65	Large					
K-81602-75	X-large					

SAFESKIN® Nitrile Exam Gloves

Purple color easily distinguishes these gloves from latex

These 100% nitrile gloves are soft, latex-free, powder-free, and nonsticky. Low chemical content reduces risk of allergic reactions. Textured fingertips enhance wet or dry grip. Beaded cuff. Ambidextrous.

Compliance: ASTM D 6319, ISO 9002.



86313-21

Specifications & Ordering Information

Catalog number	Size	Length	Thickness	Qty/bx	Price per box	
					1 bx	10 bx
K-86313-20	Small	9½"	5 mil	100		
K-86313-21	Medium			100		
K-86313-22	Large			100		
K-86313-23	X-large			90		
K-86313-30	Small	12"	5 mil	50		
K-86313-31	Medium					
K-86313-32	Large					
K-86313-33	X-large					